

# Biological Indicators & Real Time Automatic Reader

ULTRA  
RAPID  
20 MIN.

Biological Indicators (Ultra & Super rapid)  
and Incubators

Designed for Hospitals, ASC's, Physicians & Dental Practices

With increasing demands for faster sterile instruments  
turnaround in healthcare facilities, sterilization results need  
to be available as soon as possible to reliably verify the sterility  
of the load.



**Tuttnauer**  
Innovation · Legacy · Partnership

# Biological Indicators & Real Time Automatic Reader

Automatic reading and traceability system for sterilization

- 20 min or 1 hour process times
- Multiple process times can be used simultaneously
- Ease of use
- Thermal printer
- USB Connection with free software download
- One BI & reader can be used in both gravity and pre/post vacuum sterilizers



## Biological Indicators

**WTL198-0058**



**1 HOUR**

GRAVITY/VACUUM  
RAPID

**WTL198-0075**



**1 HOUR**

GRAVITY/VACUUM

**WTL198-0072**



**20 MIN**

GRAVITY/VACUUM  
ULTRA RAPID

**WTL198-0077**



**20 MIN**

GRAVITY/VACUUM


## Chemical Indicators and Tape

**WTL198-0081**



MULTI-VARIABLE INDICATOR TYPE 4

**WTL198-0082**



INTEGRATORS TYPE 5

**WTL198-0053**



SELF-ADHESIVE TAPE TYPE 1

**Tuttnauer**  
Innovation · Legacy · Partnership

Learn from our experts.

Visit our blog: [tuttnauer.com/blog](http://tuttnauer.com/blog)



Tuttnauer USA Co. Ltd.

25 Power Drive, Hauppauge, NY 11788

Tel: +800 624 5836, +631 737 4850, Fax: +631 737 0720

E-mail: [info@tuttnauerUSA.com](mailto:info@tuttnauerUSA.com), [www.tuttnauerUSA.com](http://www.tuttnauerUSA.com)

