



Clarity

DIAGNOSTICS
DIAGNOSTICS

Clarity Urocheck 120C Urine Analyzer



- ✓ Same reliability, in a small footprint analyzer
- ✓ Q.C. Lockout with Preprogrammed Values & Pass-Fail Summary
- ✓ JCAHO Compliant Hardware lockout if QC Fails
- ✓ Operator Lockout (Administrator & Operator Password Protection)
- ✓ Strip in same order as leading brand
- ✓ Fast & Easy To Use
- ✓ Patient ID Input via Manual or Bar Code Reader
- ✓ Uses Clarity Urine Strips DTG-2GP, DTG-4OB, DTG-5OB, DTG-URO7, DTG-10SG, and CD-MAC25 (Microalbumin/Creatinine)

- ✓ Large LCD Display
- ✓ Up to 120 tests per hour
- ✓ Flagging for abnormal results
- ✓ 2000 Test Internal Memory
- ✓ Automatic Paper Feeder
- ✓ Automatic Calibration and self check
- ✓ USB and RS232 Ports
- ✓ Barcode Reader (optional)
- ✓ Competitively Priced

Product Codes: CLA-U120C

Contents: 1 Urocheck 120C Urine Analyzer. Operator Manual. Quick Start Guide.

Clear Vision - Clear Results
Clarity Brand Rapid Diagnostics

For more information call 877.722.6339 or 561.288.9028
sales@claritydiagnostics.com www.claritydiagnostics.com

